
Panasonic®

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Operating Instructions Instrucciones de funcionamiento

Model No. EW3152
Modelo No. EW3152

Before operating this device, please read these instructions completely and save this manual for future use.

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Table of Contents

Introduction	3
Basics of Blood Pressure	3
Important Instructions Before Use	4
Precautions to Ensure Safe, Reliable Operation	6
Easily Check Your Blood Pressure Readings Against the JNC 7* Classification	7
Parts Identification	8
When using the AC adapter	10
When using batteries	11
Before starting measurements	12
Cautions during measurements	14
Taking a reading	14
Alarm function signals high blood pressure	17
Storing measurement values	19
Recalling stored measurement values	20
After use/Specifications	21
Cautions during measurements	22
Q & A	23
Care and Maintenance	25
Troubleshooting	26

Introduction

Thank you for purchasing the Panasonic Upper Arm Blood Pressure Monitor EW3152.

Measuring your own blood pressure is an important way of monitoring your health. High blood pressure (hypertension) is a major health problem which can be treated effectively once detected. Measuring your blood pressure between doctor visits on a regular basis in the comfort of your home, and keeping a record of the measurements, will help you monitor any significant changes in your blood pressure. Keeping an accurate record of your blood pressure will help your doctor diagnose and possibly prevent any health problems in the future.

Basic of Blood Pressure

Your heart acts like a pump, sending blood surging through your blood vessels each time it contracts. Blood pressure is the pressure exerted by blood pumped from the heart on the walls of blood vessels. Systolic pressure is the pressure exerted when the heart contracts and pumps blood into the arteries. Diastolic pressure is the pressure exerted when the heart expands, or relaxes. When you or your doctor take your blood pressure, both your systolic and diastolic pressures are measured. If your blood pressure measurement is 120 over 80 (120/80), for example, your systolic pressure is 120 while your diastolic pressure is 80.

Important Instructions Before Use

1. Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
2. If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.
3. Blood pressure can vary based on many factors, including age, gender, weight and physical condition. In general, a person's blood pressure is lower during sleep and higher when he or she is active. Blood pressure can change easily in response to physiological changes. The setting in which a person's blood pressure is measured can also affect the results. Having one's blood pressure measured by a healthcare professional in a hospital or clinic can cause nervousness and may result in a temporarily elevated reading. Because blood pressure measurements taken in a clinical setting can vary considerably from those taken at home, a person's blood pressure should be measured not only occasionally in the doctor's office, but also on a regular basis at home. Also, if you find that your blood pressure is lower at home, this is not unusual. To accurately compare with your physician's reading, take your Panasonic blood pressure monitor to your doctor's office and compare readings in this setting.
4. People suffering from cardiac arrhythmia, vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant should consult their physician before measuring their blood pressure themselves. Different values may be obtained due to their condition.
5. Try to take your blood pressure measurements at the same time and under the same conditions every day.
 - The ideal time to measure your blood pressure (to obtain your so-called "base blood pressure") is in the morning just after waking up, before having breakfast and before any major activity or exercise. If this is not possible, however, try to take measurements at a specified time prior to breakfast, and before you have become active. You should relax for about 5 minutes before taking the measurement.

- The following situations may cause substantial variations in blood pressure readings and should therefore be avoided at least 30 minutes prior to taking your blood pressure.

Blood pressure will be higher than usual:

- when you are excited or tense
- when you are taking a bath
- during exercising or soon after exercising
- when it is cold
- within one hour after eating
- after drinking coffee, tea or other beverages containing caffeine
- after smoking tobacco
- when your bladder is full

Blood pressure will be lower than usual:

- after taking a bath
- after drinking alcohol


6. Measurements may be impaired if this unit is used near a television, microwave oven, X-ray equipment or other devices with strong electrical fields. To prevent such interference, use the unit at a sufficient distance from such devices or turn the devices off.
7. This unit is designed for use by adults. Consult with your physician before using this unit on a child. Do not use on infants or toddlers.
8. This unit is not suitable for continuous monitoring during medical emergencies or operations.
9. Do not use the unit for any purpose other than measuring blood pressure. Do not use the unit together with other devices.
10. Improper handling of batteries may result in battery rupture or in corrosion from battery leakage. Please observe the following to ensure proper use of batteries.
 - a. Be sure to turn off the power after use.
 - b. Do not mix different types or sizes of batteries.
 - c. Change all batteries at the same time. Do not mix old and new batteries.
 - d. Be sure to insert batteries with correct polarity, as instructed.
 - e. Remove batteries when they are worn out, and dispose of them properly according to all applicable environmental regulations.
 - f. Do not disassemble batteries or throw them into a fire.
 - g. Do not short-circuit batteries.
 - h. Do not attempt to recharge the batteries included with the unit.

Precautions to Ensure Safe, Reliable Operation

1. Do not drop the unit. Protect it from sudden jars or shocks.
2. Do not insert foreign objects into any openings.
3. Do not attempt to disassemble the unit.
4. If the unit has been stored at temperatures below 0°C (32°F), leave it in a warm place for about 15 minutes before using it. Otherwise, the cuff may not inflate properly.
5. Do not store the unit in direct sunlight, high humidity or dust.
6. When cleaning, unplug the AC adapter from the AC outlet before cleaning. Also, do not plug in or unplug the AC adapter with wet hands. Doing so may result in electric shock or fire.
7. Do not spill water on the main body or AC adapter. Doing so may cause electric shock, short-circuiting, or product damage.
8. Do not use the AC adapter if it is damaged or if it plugs into the AC outlet loosely. Doing so may result in electric shock or short-circuiting.
9. If the product does not operate properly or seems strange, stop use and unplug the AC adapter immediately, and request inspection and service. Failure to do so could result in the risk of an accident, electric shock, or fire.
10. When using the AC adapter, make sure it is fully inserted into the AC outlet. Failure to do so could result in electric shock or short-circuiting.
11. When the product is not in use, unplug the AC adapter from the AC outlet. Failure to do so could result in deterioration of the insulation due to dust or moisture, causing a leakage electricity fire.
12. When unplugging the AC adapter, do not pull on the cord. Always grasp the AC adapter body and unplug it. Failure to do so could result in the risk of electric shock, short-circuiting, or fire.
13. Do not insert any object in the AC adapter or allow debris to collect on it. Doing so could result in the risk of electric shock, short-circuiting, or fire.
14. Do not scratch, damage, modify, forcibly bend, pull, or twist the AC adapter cord. In addition, do not place heavy objects on top of the cord or pinch it with objects. Do not wrap the cord around the main body and the AC adapter. Doing any of these things may result in the risk of fire or electric shock.
15. Do not step on the AC adapter. Doing so may cause an accident or injury.
16. Do not use a voltage converter. Doing so may cause fire or electric shock.
17. This AC adapter is designed for use only with Panasonic Blood Pressure Monitors, model EW3152. Use with other electrical devices may result in death or serious personal injury.
18. ALWAYS unplug the AC adapter from the power outlet before performing maintenance.
19. Do not let this product to be used by children.
20. Do not let this product to be used by people who cannot express their feelings, or who cannot operate the product by themselves. Such use may result in an accident or injury.
21. Use AC adapter only with electrical outlet of voltages marked on name plate. For use in the U.S.A., use AC adapter only with AC 120V electrical outlet. Use with electrical outlets at other voltages may results in the risk of fire, electric shock, or burns.
22. Do not allow any material to drape over the AC adapter. Doing so could result in the risk of fire.
23. Do not modify or attempt to repair the product or AC adapter by yourself. Doing so could result in the risk of fire, electric shock, or personal injury.

Easily Check Your Blood Pressure Readings Against the JNC 7* Classification

Blood Pressure Categories

If your systolic pressure falls in one category but your diastolic pressure in another, your level is classified in the higher of the two categories.



* JNC 7: The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation and Treatment of High Blood Pressure National Institute of Health (NIH) Publication; No. 03-5233, May 2003


IMPORTANT:

- Do not be alarmed by temporarily high or low readings because fluctuations in a person's blood pressure are not uncommon. If possible, measure and record your blood pressure at the same time every day, and consult your physician if you have questions or concerns.
- If abnormal variations in blood pressure are observed in measurement, please consult your physician.



PARTS IDENTIFICATION

Please check the instrument.


Back side of main body


Included accessory

AC adapter

WHEN USING THE AC ADAPTER

* Use only the exclusive AC adapter (included accessory). (See p. 9.)

- 1 Insert** the AC adapter plug into the connector on the side of the main unit.
- 2 Plug** the AC adapter fully into the AC outlet. (An operating sound may be heard momentarily when the AC adapter is plugged into the outlet.)

Be sure to use only the exclusive AC adapter (included accessory).

Use of other AC adapters could cause an accident.


Caution when using only the AC adapter

When no batteries are installed in the main unit and the AC adapter is unplugged from the AC outlet or main unit, the set date and time will be reset and deleted. Because of this, even when using the AC adapter, it is recommended that batteries be installed in the main unit.

When the AC adapter is used, power is supplied from the AC adapter, so virtually no power from the batteries is used.

WHEN USING BATTERIES

* Always use alkaline batteries. (Four AA-size LR6 alkaline batteries)
(Use only alkaline batteries. Do not use other types of batteries.)

1 Turn the main body over, press lightly on the \triangle mark on the battery cover, and **remove** the cover by sliding it in the direction of the arrow.

2 **Slide** the battery holder **out** in the direction of the arrow.

3 Check the polarities of the batteries and **install** them correctly.

4 **Slide** the battery holder back **in**. (An operating sound may be heard momentarily when the battery holder is put back in.)

5 Replace the battery cover and slide it **closed**.

Replace batteries when:

- The mark shown at left appears during use.
- All displays are lit for a moment after the On/Off Start button is pressed, and then the display goes blank.
- No display appears even if the On/Off Start button is pressed.

- Approximately 200 measurements (approximately 2 months when measurements are taken 3 times per day) can be taken with a set of fresh Panasonic alkaline batteries (four AA-size LR6 batteries). (Measurement conditions: Room temperature 22°C (72°F); Pressurized to 170mmHg; Upper arm circumference: 30cm (11-51/64")) (The supplied batteries are included for trial purpose only. Accordingly, there is no guarantee regarding the length of their useful lives.)
- Even if the product is not used for a long time, be sure to replace batteries once a year. If battery replacement is delayed, damage due to battery leakage may occur.
- When replacing batteries, always use four fresh alkaline batteries of the same type from the same manufacturer.
- When room temperature is low, battery life may become shorter.

BEFORE STARTING MEASUREMENTS

Open the display

Press the Open button for display on the front of the main unit.


* The display will only open partially by itself. Use your hand to open the display fully for use.

Setting the date and time

* Please be sure to set the date and time before taking the first measurement after purchasing the blood pressure monitor or after changing the batteries.

- If measurements are taken without setting the date and time, the date and time will not be stored.
- Please read the "Caution when using only the AC adapter" on p. 10.

1 Press the Set button below the display. The month value will blink.

8/ 3 0:00

4 Set the hour.

- ① Press the **Set button**. The hour value will blink.
- ② Press the **Adjust button** to set the correct hour.

8/ 1 0:00

2 Set the month. Press the **Adjust button** repeatedly to set the correct month.

8/ 3 15:00

5 Set the minute.

- ① Press the **Set button**. The minute value will blink.
- ② Press the **Adjust button** repeatedly to set the correct minute.

8/ 1 0:00

3 Set the day.

- ① Press the **Set button**. The day value will blink.
- ② Press the **Adjust button** repeatedly to set the correct day.

8/ 3 15:30

6 Press the Set button at the exact time the minute changes **to finalize** the date/time setting.

7 Press the On/Off Start button to switch off the power.

When the date or time are blinking, the clock is stopped.

During leap years, 3/1 will be displayed on February 29. For this reason, reset the date to 2/29.

Use the tip of a ballpoint pen or similar object to lightly press the Set and Adjust buttons. If you use a mechanical pencil, the lead may break, causing a malfunction.

Getting into the measurement position

Either the left or right arm can be measured. However, when measurements are taken using the right arm, it may be difficult to see the display.

* Since there may be blood pressure differences of around 10mmHg between the left and right arm, be sure to use the same arm for all measurements.

* The product cannot be used while lying down. Please sit upright during use.

Place the main unit on a table or similar object so that the main unit is at the same height as your heart.

* Differences in height will cause differences in blood pressure values.

If the heights of the chair and table are not suitable, you may have to lean forward. In this case, form a stable stack of books or similar objects under the main unit to adjust the height of the main unit.

1 Sit next to the table as shown in the figure at left.

2 Slide your arm through the measuring section.

Pull the main unit close to you so that you do not have to lean to reach it.

- Remove your shirt or roll up your sleeve so that the measuring section will be on your bare skin.
- If rolling up your sleeve causes pressure on your arm or underarm, remove your shirt for measurements.



Be sure that the rolled-up sleeve does not get inside the measuring section when taking measurements.

3 Place your elbow **in the support.**

- If it is difficult to find the elbow support, use your other hand to locate it.

4 Stretch out your arm and **sit straight.**

- The unit will naturally pivot upwards up to 15° as you place your hand on the table in a relaxed position.

CAUTIONS DURING MEASUREMENTS

In the following situations, accurate measurements cannot be taken.

- Within 1 hour of eating or drinking alcohol
- Immediately after drinking coffee or tea, or smoking
- Immediately after exercising or taking a bath
 - Wait at least 20 minutes and take a measurement while staying quiet.
- In a moving vehicle
- In cold places
 - Take measurements at ambient temperatures of 10°C to 40°C (50°F to 104°F).
- When you need to urinate
 - Wait several minutes after urinating before taking measurements.

TAKING A READING

Cautions to obtain accurate measurements
Sit quietly and relax while taking measurements.

Do not lean forward with your body at an angle during measurements.

Take measurements in the proper posture.

Do not take measurements at the elbow. Place your elbow in the elbow support.


Do not lift the main unit off the table. Keep it on the table during measurements.

Do not lift your arm. Rest your hand on the table with your palm upwards.

Do not let your rolled-up sleeve be inside the measuring section. Roll your sleeve far up your arm.

Be sure to obey the following points for accurate measurements:

- Do not move your body, arm, or fingertips.
- Do not touch the main unit.
- Do not talk.
- Do not use a cellular phone near the main unit.


Not following these points may make it impossible to obtain accurate measurements.

The angle of this product can be adjusted so that measurements can be taken without having to be in a forced posture.

- The angle can be adjusted only when the display is open.
- The angle is not locked.

Even if you sit still, the angle may change slightly during pressurization.

However, there is no effect on accuracy even if the angle changes during measurements.


TAKING A READING (CONTINUED)

Measurement

1 Press the On/Off Start button.

- ① The pressure cuff will be tightened automatically.
- ② Pressurization will be performed automatically and measurements will be taken.

The entire display will light for approximately 1 second.

The pressure cuff will be tightened automatically.



The heart mark ♥ will blink when pulse is detected during pressurization.

2 When measurement has been completed, **the blood pressure and pulse values will be displayed.**

After measurement values have been shown for approximately 5 seconds, the pressure cuff will be loose. (However, the time depends on the size of the arm.)

- The blood pressure indicator lamps will light (for approximately 12 seconds).
- If the measured value is in the high blood pressure region, the values will blink (for approximately 12 seconds). (For how to read the blood pressure display, see p. 17.)
- The ① ② marks will light. (For storing measurement values, see p. 19.)
- * The pulse rate value is the beats per minute calculated based on the pulse rate during measurement.

3 Remove your arm from the main unit.

4 To switch off the blood pressure monitor, **press** the On/Off Start button.


- If you forget to switch it off, the blood pressure monitor will switch off automatically after approximately 5 minutes.

Alarm function signals high blood pressure


Blood pressure values for both systolic and diastolic readings are determined and indicated via a three-color LED system as being within normal or hypertensive ranges based on definitions and classification of blood pressure levels by the JNC 7.

Moreover, when blood pressure is in the high range, the reading displayed on screen flashes on and off to alert user.

According to JNC 7 Classification, values consistently in excess of 140 mmHg (Sys.) and/or 90 mmHg (Dia.) are considered to constitute high blood pressure.

- If the measured blood pressure value is in the high blood pressure region, the values will blink for approximately 12 seconds. (Only immediately after measurement)

- **Systolic blood pressure:** 140mmHg or higher
- **Diastolic blood pressure:** 90mmHg or higher

If measurements could not be taken properly

The remeasure mark E will be displayed.

When taking successive measurements, press the On/Off Start button to switch off the blood pressure monitor once and stay quiet for 4 to 5 minutes, and then take the next measurement.

Advice:


In general, blood pressure measurements are lower for the second measurement because you are less nervous and are more tranquil.

For people who easily get tense, measuring twice each time is very effective.

When two measurements are taken, normally the smaller measurement is stored.

When taking two measurements, always stay quiet for 4 to 5 minutes after the first measurement before taking the second measurement.


To interrupt measurements quickly


If the pressure cuff will not deflate even if the On/Off Start button is pressed, press the Quick air release button and hold it pressed. The air will be released in approximately 4 seconds.

If the pressure cuff comes off:

- 1** With the belt and seam downwards, **slide** the pressure cuff into the main body.


- 2** **Press** the pressure cuff **into** the groove of the main body. Either the front or the back can be attached first.

* When re-attaching, be sure not to touch the inside of the main body where the pressure cuff has been removed. Doing so may result in an accident or injury.

STORING MEASUREMENT VALUES

This product can separately store 90 measurements each for two persons.

After measurement, the ① ② will light.


Press the M/R (memory/recall) button with the number corresponding to the person for which measurements should be stored.

Example: Person A: Press ①

Person B: Press ②

Storing is completed.

- Up to 90 measurement values can be stored. When the memory is full and additional values are stored, the oldest value in memory will be deleted and the new value will be stored.
- If measurements could not be taken properly (the E mark is shown instead of the blood pressure value), measurement values cannot be stored.
 - * However, if the E mark is displayed for only the pulse value, data can be stored.

Number of stored data

To erase all measured values stored in memory:

1. Recall the measurement value for the person whose data you want to delete by pressing that person's M/R (memory/recall) button.
2. Press the same M/R (memory/recall) button again, and hold it pressed until the measurement values disappear (approximately 3 seconds or more).

* It is not possible to erase single measurement values.


To switch off the blood pressure monitor,

press the On/Off Start button.

- If you forget to switch it off, the blood pressure monitor will switch off automatically after approximately 5 minutes.


RECALLING STORED MEASUREMENT VALUES

- 1 Press** the M/R (memory/recall) button for the person whose data you want to recall.

The average of all measured values in memory will be **shown**.



- **Avg.** mark will be shown at the bottom of the display.
- For example, if 90 measured values are stored, the average of the 90 measurements will be shown. (If only 1 measurement is stored, the **Avg.** mark will not be shown.)
- The blood pressure indicator lamps will light (for approximately 6 seconds).

- 2 Press** the M/R (memory/recall) button again.

Each time you press the M/R (memory/recall) button, the previously stored measurement values will be **shown** in order, starting with the newest data.


- The blood pressure indicator lamps will light (for approximately 6 seconds).
- If no measurement values are stored, no value will be shown.


Press the M/R (memory/recall) button lightly for only a short time. If you hold the button pressed for approximately 3 seconds or longer, all data in memory will be erased, so please be careful.

- 3 To cancel measurement value recall,**
press the On/Off Start button.

- If you forget to switch it off, the blood pressure monitor will switch off automatically after approximately 30 seconds.

AFTER USE/SPECIFICATIONS

1 If the main body is set at an angle, **return it to normal.**

2 Fold the display up against the main body for **storage.**

- Fold it up until you can feel a click with your hand.

Specifications (Main body)

Power source	DC 6V (4 LR6 AA-size alkaline batteries) AC 100-240V 50-60Hz (When using AC adapter)
Display	Digital LCD
Measurement method	Oscillometric system
Measurement range	Blood pressure: 0 to 280mmHg Pulse rate: 30 to 160 beats per minute
Accuracy	Blood pressure: Within ± 3 mmHg Pulse rate: Within $\pm 5\%$
Operation temperature / humidity range	10 to 40°C (50 to 104°F), 15 to 90% RH
Storage temperature/ humidity range	-20 to 60°C (-4 to 140°F), 10 to 95% RH
Measurable arm circumference	Approx. 20 to 34cm (7-3/4" to 13-1/4")
Weight	Approx. 4.9 lbs. (1200g) (not including batteries)
Dimensions	22 × 15.5 × 21cm (8-21/32" × 6-3/32" × 8-1/4")


Specifications (AC adapter)

Input	AC 100-240V, 50-60Hz
Output	DC 6V 2A

CAUTIONS DURING MEASUREMENTS

In the following situations, accurate measurements cannot be taken.

Within 1 hour of eating or drinking alcohol

Immediately after drinking coffee or tea, or smoking

Immediately after exercising or taking a bath

In a moving vehicle

Wait at least 20 minutes and take a measurement while staying quiet.

In cold places

Take measurements at ambient temperatures of 10 to 40°C (50 to 104°F).

When you need to urinate

Wait several minutes after urinating before taking measurements.

Q & A

Q Why are the blood pressure values measured at home different from those measured at the doctor's office?

A Blood pressure varies 24 hours a day.

In addition, it can change greatly depending on the weather, your mental state, exercise, etc.

In particular, there is the phenomenon of "White uniform hypertension" at doctor's offices, where the blood pressure measured in a doctor's office is higher because of stress or worry.

In addition, if the blood pressure measured at home is sometimes higher and sometimes lower than the blood pressure measured at the doctor's office, the following can be considered. Please pay attention to these points when taking measurements.

① **Measuring posture was poor.**

Did you take measurements while leaning forward with your body at an angle?

② **Measurements were not taken at the upper arm.**

Please take measurements with your elbow in the elbow support of the main body. (See p. 13.)


③ **Were you worried or upset during measurements?**

Before taking measurements, take 2 or 3 deep breaths and relax, so that your blood pressure stabilizes. In addition, it is recommended that you sit quietly for about 5 minutes before taking measurements.

Q Why are measured blood pressure values different each time I take a measurement?

A

- ① Blood pressure changes even within a single day, and it is different depending on the measurement posture. Because of this, when taking measurements, be sure to always take measurements under the same conditions.
- ② For persons taking blood pressure medication, blood pressure may change greatly due to the effects of the medication.
- ③ When taking 2 measurements in succession, be sure to stay quiet for 4 to 5 minutes after the first measurement before taking the next measurement.

Q Why are the measured values different for the left and right arm?

A

It is said that differences in blood pressure values between the left and right arms may be around 10mmHg. The blood pressure can be measured in either the left or right arm, but all measurements should be taken using the same arm for comparative consistency.

CARE AND MAINTENANCE

Do not apply excessive force to the instrument or drop the instrument.

- Doing so may cause damage.

If the instrument is stored at temperatures below the freezing point, do not use it immediately.

Leave it in a warm place for at least an hour before use.

- If the instrument is not allowed to warm up, it may not pressurize.

Do not insert dust or foreign objects into the instrument.

- Doing so may result in damage.

If the instrument becomes dirty, clean it with a soft cloth moistened with warm water or soapy water.

(Do not use alcohol, benzene, or paint thinner.)

- Use of such chemicals may result in cracking or discoloration.

For storage, avoid high temperatures, high humidity, and direct sunlight.

- Storage under such conditions may result in damage.

The pressure cuff cannot be washed.

TROUBLESHOOTING

Display	Status before error	Cause and solution	If problem continues
<p>Remeasure mark E is displayed.</p> 	<p>Pressurization to more than 280mmHg was performed.</p> <p>Pressure decreased suddenly and the remeasure mark E appeared.</p> <p>The heart mark ♥ blinked only a few times or not at all.</p> <p>No pressurization</p>	<p>Were measurements taken according to the proper procedure and in the correct posture? See p. 13 and 14.</p> <p>Be sure your arm is properly inserted into the main unit. See p. 13.</p>	
Symptom	Possible cause		
Systolic or diastolic blood pressure value is high or low.	<ul style="list-style-type: none"> You moved or talked during measurements. Your arm was not inserted properly into the main unit. <p>See p. 13 and 14.</p>	<p>If problem continues</p> <p>Ask your doctor for advice.</p>	
Value is different from that measured by the doctor.	<ul style="list-style-type: none"> Blood pressure varies in response to minute changes in your mental state, such as your reaction to a doctor's presence. Relax for 4 to 5 minutes and take another measurement. <p>See p. 14 and 17.</p>		
Measurement values are different with each measurement.	<ul style="list-style-type: none"> Batteries are installed with the polarities wrong. <p>See p. 11.</p>		
No display appears even when the On/Off Start button is pressed.	<ul style="list-style-type: none"> Batteries are exhausted. <p>See p. 11.</p>		
When the On/Off Start button is pressed, the display disappears during pressurization (⊕⊗ mark appears).	<ul style="list-style-type: none"> Since measurements are taken during pressurization, operating sound may be heard and pressurization rate may change during pressurization. This is normal and is not a malfunction. 		
During pressurization, the sound or pressurization rate changes.	<ul style="list-style-type: none"> Press the ring portion on the outside of the pressure cuff into the main body. <p>See p. 18.</p>		
Pressure cuff comes off.			
Instrument is broken.	<p>Take it to the store where it was purchased for testing and repairs.</p>		

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Division of Panasonic Corporation of North America**